

# Technical note

## BAKERBOND™ Multimode ion exchange chromatography media

### BAKERBOND™ POLYABX

#### INTRODUCTION

PolyABx is a multimode ion exchanger which primarily functions as a weak cation exchanger over a wide pH range. Secondary anion exchange sites are due to the presence of amine groups. This multimodal functionality offers better selectivity than conventional weak cation exchange media with equivalent capacities and is capable of separating proteins and peptides having similar isoelectric points (pI). Unique selectivity of closely related molecules is often achieved with PolyABx where conventional ion exchangers fail to provide sufficient separation in a process environment.

The average particle size of 35 µm with narrow particle size distribution ( $d_{60}/d_{10} < 1.4$ ) provides high efficiency and will produce narrower elution bands of highly concentrated product. This increases resolution and decreases pool volumes compared to media based on larger particles. High efficiency in combination with high selectivity enables higher loading levels while maintaining separation.

PolyABx is based on highly cross linked rigid spherical polymethacrylate particles with optimal porosity and mechanical strength. The hydrophilic backbone has low non-specific binding and the ability to withstand prolonged contact with all commonly used cleaning and sanitizing solutions. PolyABx can be easily packed to bed heights up to 40 cm and operated in the linear velocity range of up to 500 cm/hr using conventional columns with pressure ratings of 2 to 10 bar. PolyABx is manufactured under the strictest controls and testing and supported by world class quality systems and applications resources.


As a result of this unique combination of selectivity, efficiency, and the ability to maintain performance at high linear velocities, PolyABx can provide high process throughput, productivity improvements, and enable new processes.

#### Characteristics

Functionality	Primarily a weak cation exchanger with weak anion exchange sites	
Functional Group	-NH-C(=O)-CH <sub>2</sub> CH <sub>2</sub> COOH	
Ion Exchange Capacity	Cation exchange: 0.15 – 0.25 H+meq/mL	Anion exchange: 0.10 – 0.20 Cl meq/mL
Median Particle Size	35 µm	
Particle Size Distribution ( $d_{60}/d_{10}$ )	< 1.4	
Median Pore Size	500 Å	
Exclusion Limit	1 × 10 <sup>6</sup> Daltons	
Operating pH Range	4.5 – 14.0	
Cleaning pH Range	1.0 – 14.0	
Chemical Stability	All commonly used aqueous buffers, sodium hydroxide, acetic acid, phosphoric acid, guanidine hydrochloride, up to 100% ethanol, methanol, or 2-Propanol.	
Shipping Solvent	Media is shipped as 1:1 slurry in acetate buffer containing benzyl alcohol at pH 4.5 and can be stored for up to 5 years at 4 – 15 °C. The media can be stored in 20% ethanol or 0.1 M NaOH for up to 5 years.	

## MULTIMODE FUNCTIONALITY


Fig. 1. Ligand structure of PolyABx



The multimode functionality of PolyABx is obtained by covalently bonding PEI to the surface of highly cross linked polymethacrylate beads and functionalizing PEI to obtain carboxylic acid groups, thereby providing a primarily weak cation exchange and secondary weak anion exchange behavior.


## PHYSICAL CHARACTERISTICS

Fig. 2. SEM image of PolyABx at 156 magnification




Uniform particle size distribution is evidenced in this picture.

Fig.3. SEM image of PolyABx at 7600 magnification



The picture shows excellent and uniform pore structure of the bead.

Fig.4. Particle size distribution of a typical lot of polymethacrylate beads




The median particle size of 35 micron and a particle size distribution of less than 1.4 ( $d_{90}/d_{10}$ ) shows a total lack of fines.

As illustrated by the SEM picture and particle size distribution graph, PolyABx provides high column efficiency leading to narrow elution bands and therefore concentrated product fractions. The SEM picture also shows good porosity and has an average pore size of 500 Å for an exclusion limit of  $1 \times 10^5$  dalton. This open porosity coupled with the bead's high mechanical strength allows the use of PolyABx at high linear velocities due to good mass transfer.

## SEPARATION PERFORMANCE: SELECTIVITY AND CAPACITY

**Fig. 5.** Separation of bovine serum albumin (bsa), igg, cytochrome c and lysozyme



The unique selectivity of PolyABx is capable of separating proteins with closely related isoelectric points.

**Sample injection:** 100 mL

**Sample:** 1B SA( pI4 .7), 2I gG( pI7 .2), 3C ytochrome( pI1 0.4), 4L ysozyme( pI1 1)were dissolved in binding buffer

**Column:** Millipore V90 165 × 90 mm ID packed with PolyABx


**Binding Buffer:** 20 mM Sodium Acetate/Acetic Acid, pH 5.5.

**Elution Buffer:** 1.0 M NaCl in binding buffer

**Flow rate:** 300 cm/hour.

**Gradient:** 10 CV

**FIG. 6.** Dynamic binding capacity of rabbit igg on polyabx at different linear velocities



Dynamic binding capacity is maintained at high linear velocities.

**Rabbit IgG concentration:** 1 mg/mL in Buffer A

**Column:** 4.6 × 100 mm.

**Binding buffer:** 20 mM sodium acetate at pH 5.6.


**Elution buffer:** 1 M acetate at pH 6.0

PolyABx media provides unique selectivity and can be used to separate closely related proteins due to its unique chemistry and multimodal functionalities. As Fig. 5. illustrates, Lysozyme and Cytochrome having close isoelectric points (pI 10.4 and 11) are well resolved and can be easily separated. As expected, BSA flows through without binding while IgG is retained demonstrating the capability to separate albumin and IgG as well. This in combination with the media's high

efficiency enable high column loading. Fig. 6. Illustrates that the breakthrough capacity of IgG does not change significantly over the linear velocity range of 180 to 720 cm/hour. This is due to the particle's open pore structure and rigidity facilitating good mass transfer. High resolution and the ability to maintain capacity at high linear velocities facilitate high throughput for both large and small proteins.


## STABILITY AND CONSISTENCY

Fig. 7. Polyabx performance consistency over 50 cycles




Separation performance is maintained over 50 cycles.

Fig. 8. Dynamic binding capacity of rabbit igg on polyabx under various storage conditions



Dynamic binding capacity is maintained under various storage conditions.

Fig. 9. Separation of Bovine Serum Albumen, IgG and Lysozyme.



Separation performance is maintained under various storage conditions.

In Fig. 7., the consistent separation behavior of PolyABx was evaluated by separating an IgG and Lysozyme mixture for 50 cycles. After each separation the column was regenerated by washing with 5 CV of elution buffer. There is no measurable change in retention time of the two proteins, indicating consistent separation performance over multiple cycles

**Sample:** 0.25mg IgG and 0.12mg lysozyme in binding buffer;  
**Sample injection:** 0.5mL

**Column:** 10 × 10 mL packed with PolyABx

**Conditions:** Linear gradient from binding buffer (0.05M MES/NaOH pH 5.6) to elution buffer (1 M NaCl in 0.05 M MES/NaOH pH 5.6) in 5 column volumes (26 min)

**Flow rate:** 1.5 mL/min

**Blue:** 1st injection

**Red:** 20th injection

**Pink:** 50th injection

**Rabbit IgG concentration:** 1 mg/mL in buffer A

**Column:** 4.6 × 100 mm.

**Binding buffer:** 20 mM sodium acetate at pH 5.6.

**Elution buffer:** 1 M acetate at pH 6.0

**A:** Initial

**B:** 24 months in storage buffer at 4 °C

**C:** 5 years in 0.1 N NaOH at 4 °C

**D:** 5 years in 0.1 N NaOH at room temperature

**E:** 5 years in storage buffer at 4 °C after washing with 0.5N NaOH

**Column:** 0.775 × 10 cm packed with PolyABx

**Buffer A:** 50 mM MES, pH 6.2 and Buffer B: 50 mM MES containing 1M NaCl, pH 6.2.

**Flow rate:** 4 mL/min

**Orange:** 20% ethanol at room temperature, Pink: Phosphoric acid at room temperature


**Yellow:** Phosphoric acid at 4 °C, Green: Sodium hydroxide at room temperature

**Blue:** Sodium hydroxide at 4 °C

PolyABx is chemically stable and as shown in Figs 8 and 9., there is no significant change in capacity or resolution when media is subjected to various storage conditions. The former study evaluated the change in dynamic binding capacity of human IgG between initial media compared to media samples stored in different solutions for up to 5 years. Slight variations are expected within the error of the experiment. High selectivity and resolution of PolyABx are also maintained when stored at extreme pH conditions over a period of 10 months as illustrated in the 5 overlaid chromatograms where 5 different storage buffers were chosen to demonstrate consistent chromatographic behavior.

## OPERATIONAL FLEXIBILITY

Fig. 10. Efficiency of column (20 cm id x 18 cm high) packed with PolyABx




13,926 plates/meter is achieved with an asymmetry of 1.18

Fig. 11. Efficiency of column (45 cm id x 18 cm high) packed with PolyABx



13,500 plates/meter is achieved with an asymmetry of 1.20.

Fig. 12. Pressure–Flow relation



The linear behavior of the pressure-flow curve indicates that there is no resin compression.

Repeated column packing of the media yields high efficiencies as well as low asymmetries and is reproducible from a 20 cm id to 45 cm id column, indicating easy process scalability.

The narrow particle size distribution and mechanical stability of PolyABx provides the ability to operate in conventional low pressure columns at a bed heights of up to 20 cm. Modern medium pressure columns, rated at 7 bar, can be used to pack beds up to 40 cm bed height and operate at higher linear velocities. The linear pressure–flow curve for PolyABx indicates that there is no resin compression at a bed height of 18 cm and a column diameter of 20 cm. In addition, the media does not need to be defined prior to the initial packing or subsequent packings, reducing the consumption of buffers and simplifying operations.

Fig. 13. Scale-up to a BPG column



PolyPEI can easily be packed in large columns for high throughput applications.

### Product ordering information

Size	Product Number
50 mL	7586-01
100 mL	7586-02
500 mL	7586-03
1 L	7586-04
5 L	7586-05

As our channel brand, VWR, part of Avantor, offers an integrated, seamless purchasing experience that is optimized for the way you do business so you can obtain what you need to move your processes forward. For more information or to buy our products, please visit your local VWR website.

#### AUSTRIA

VWR International GmbH  
Graumannsgasse 7  
1150 Vienna  
Tel.: +43 1 97 002 0  
info.at@vwr.com

#### BELGIUM

VWR International bv  
Researchpark Haasrode 2020  
Geldenaaksebaan 464  
3001 Leuven  
Tel.: +32 (0) 16 385 011  
vwr.be@vwr.com

#### CANADA

VWR International  
2360 Argentia Road  
Mississauga, Ontario L5N 5Z7  
Tel.: +1 800 932 5000  
Canada\_Orders@vwr.com

#### CHINA

VWR (Shanghai) Co., Ltd  
Bld.No.1, No.3728 Jinke Rd,  
Pudong New District  
Shanghai, 201203- China  
Tel.: 400 821 8006  
info\_china@vwr.com

#### CZECH REPUBLIC

VWR International s. r. o.  
Veetee Business Park  
Pražská 442  
CZ - 281 67 Stříbrná Skalice  
Tel.: +420 321 570 321  
info.cz@vwr.com

#### DENMARK

VWR International A/S  
Tobaksvejen 21  
2860 Søborg  
Tel.: +45 43 86 87 88  
info.dk@vwr.com

#### FINLAND

VWR International Oy  
Valimotie 9  
00380 Helsinki  
Tel.: +358 (0) 9 80 45 51  
info.fi@vwr.com

#### FRANCE

VWR International S.A.S.  
Le Périgares – Bâtiment B  
201, rue Carnot  
94126 Fontenay-sous-Bois cedex  
Tel.: 0 825 02 30 30\* (national)  
Tel.: +33 (0) 1 45 14 85 00  
(international)  
info.fr@vwr.com  
\* 0,18 € TTC/min + prix appel

#### GERMANY

VWR International GmbH  
Hilpertstraße 20a  
D - 64295 Darmstadt  
Tel.: 0800 702 00 07\* (national)  
Tel.: +49 (0) 6151 3972 0 (international)  
info.de@vwr.com  
\*Freecall

#### HUNGARY

VWR International Kft.  
Simon László u. 4.  
4034 Debrecen  
Tel.: +36 52 521130  
info.hu@vwr.com

#### INDIA

VWR Lab Products Private Limited  
No.139, BDA Industrial Suburb,  
6th Main, Tumkur Road, Peenya Post,  
Bangalore, India – 560058  
Tel.: +91 80 28078400  
vwr\_india@vwr.com

#### IRELAND

VWR International Ltd  
Orion Business Campus  
Northwest Business Park  
Ballycoolin  
Dublin 15  
Tel.: +353 (0) 1 88 22 222  
sales.ie@vwr.com

#### ITALY

VWR International S.r.l.  
Via San Giusto 85  
20153 Milano (MI)  
Tel.: +39 02 3320311  
info.it@vwr.com

#### KOREA

VWR International ~  
17, Daehak 4-ro, Yeongtong-gu,  
Suwon-si, Gyeonggi-do  
Tel.: +82 31 645 7256  
saleskorea@avantorsciences.com

#### THE NETHERLANDS

VWR International B.V.  
Postbus 8198  
1005 AD Amsterdam  
Tel.: +31 (0) 20 4808 400  
info.nl@vwr.com

#### MEXICO

VWR International, S.de R.L.  
de C.V.  
Km. 14.5 Carretera  
Tlalnepantla-Cuautitlán  
Col. Lechería  
Tultitlán Edo. de México  
CP 54940  
Tel.: +52 (55) 5005 0100  
vwrmx@vwr.com

#### MIDDLE EAST & AFRICA

VWR International FZ-LLC  
Office 203, DSP Lab Complex,  
Dubai Science Park,  
Dubai, United Arab Emirates  
Tel.: +971 4 5573271  
Info.mea@vwr.com

#### NORWAY

VWR International AS  
Brynsalleen 4,  
0667 Oslo  
Tel.: +47 22 90 00 00  
info.no@vwr.com

#### POLAND

VWR International Sp. z o.o.  
Limbowa 5  
80-175 Gdansk  
Tel.: +48 58 32 38 200  
info.pl@vwr.com

#### PORTUGAL

VWR International -  
Material de Laboratório, Lda  
Centro Empresarial de Alfragide  
Rua da Indústria, nº 6  
2610-088 Amadora  
Tel.: +351 21 3600 770  
info.pt@vwr.com

#### SINGAPORE

VWR Singapore Pte Ltd  
18 Gul Drive  
Singapore 629468  
Tel.: +65 6505 0760  
sales.sg@vwr.com

#### SPAIN

VWR International Eurolab S.L.U.  
C/ Tecnología 5-17  
A-7 Llinars Park  
08450 - Llinars del Vallès  
Barcelona  
Tel.: +34 902 222 897  
info.es@vwr.com

#### SWEDEN

VWR International AB  
Fagerstagatan 18a  
163 94 Stockholm  
Tel.: +46 (0) 8 621 34 00  
kundservice.se@vwr.com

#### SWITZERLAND

VWR International GmbH  
Lerzenstrasse 16/18  
8953 Dietikon  
Tel.: +41 (0) 44 745 13 13  
info.ch@vwr.com

#### UK

VWR International Ltd  
Customer Service Centre  
Hunter Boulevard - Magna Park  
Lutterworth  
Leicestershire  
LE17 4XN  
Tel.: +44 (0) 800 22 33 44  
uksales@vwr.com

#### UNITED STATES

VWR International, LLC  
100 Matsonford Road  
Building One Suite 200  
Radnor, PA 19087  
Tel.: +1 800 932 5000  
VWRCustomerService@vwr.com